

eFlexx

Compact Power Mobility

Adjustable, compact, lightweight, and stylish, eFlexx is the perfect general-use power wheelchair to extend your journey and take you anywhere!



Self-Propel

20" rear wheels with freewheel hubs allow users to easily switch between power drive and self-propelled modes.

Quick-Release Battery Pack

The 26 Ah lithium battery provides up to 20 km* of driving range.



Transport

14" rear wheels offer the lightest weight and smallest storage size



Compact for Travel

In addition to having a foldable frame, eFlexx's battery pack and footrests are removable, making it compact enough to fit into a car.

Flexible and Adjustable Design

Configurable seat width and depth, swing-away footrests, backrest angle, and armrest height are all adjustable for a personalized fit.



Benefits of Lithium Batteries

Lightweight, environmentally-friendly, fast-charging, and long-lasting.



Swing-away footrests

Enable forward transfers into and out of the wheelchair



Flip-up armrests

Allow users to laterally transfer into and out of the wheelchair.



Height adjustable armrests

Tool-free height adjustment from 210 - 290 mm to suit user needs.

Foldable push handles

Allow attendants to push the wheelchair and easily fold away when not in use.



Swing-away controller bracket

Approach tables with ease, without damaging your controller or injuring your hand.

User-friendly design

The intuitive, high visibility adjustment levers allow users to easily adjust the armrests and footrests.

OPTIONS

- Wheel Locks
- Rigidifying push bar
- Vertical footrests
- Elevating legrests
- Swing-away controller bracket
- Foldable headrest

ACCESSORIES

- Handle bar
- Cane holder
- IV pole

Product Name	Nominal seat width (mm)	Nominal seat depth (mm)	Effective seat width (mm)	Effective seat depth (mm)	Backrest height w/Cushion (mm)	Seat height w/o Cushion (mm)	Seat height w/Cushion (mm)	Power seat lift (mm)	Power Tilt, Seat Plane (X°)	Backrest recline, seat to back angle (X°)	Legrest angle, leg to seat (X°)	F/M/R Wheel (in)	Overall length with legrest (mm)
eFlexx-14	360/410/460	415/440/465	405/455/500	560/585/610	410	-	485	-	-	90-110	99	8/-/14	1015-1065
eFlexx-20	360/410/460	415/440/465	390/455/510	545/585/610	410	-	485	-	-	90-110	99	8/-/20	1125-1175

Overall width (mm)	Max. User weight (kg)	Drive range (km)	Max. speed (km)	Obstacle climbing (mm)	Min. turning radius w/o user (mm)	Ability to climb rated slope (X°)	Battery (Ah)	Electrical system
560/615/665	100	20	8	35	795/825/905	10	26.1*	Linx
590/645/690	100	20	8	35	795/840/865	10	26.1*	Linx

- The actual product specifications may differ from the data in the table by ±10 mm ±0.5 kg.
 - The actual product specifications may vary according to different configurations.
 - Karma reserves the right to modify information herein without further notice.

