



PATENTED

# Ergo Anti-Tipper

Overcome Curbs With a Single Step




## Curb-Climbing Mode


Step on the tipper pedal to automatically raise the anti-tipper out of the way for curb climbing.

## Safety Mode

Step on the tipper head to lower and automatically lock the anti-tipper in the safe operating position.



Before climbing curb, step on the tipping lever to quick-release the anti-tipper.



After climbing a curb, make sure you push down the tipper to prevent wheelchair from tipping over.

### Standard For Models

Flexx Adapt

### Optional For Models

Ergo 2  
Ergo 3  
Ergo Lite  
Ergo Lite 2



Scan to learn more about Ergo Anti-tipper



Why do you need the Ergo Anti-tipper?



How to Install the Ergo Anti-tipper  
Rev.01 Mar. 2025

KARMA MEDICAL PRODUCTS CO., LTD.

T +886-5-206 6688

F +886-5-206 7788

Email : [globalsales@karma.com.tw](mailto:globalsales@karma.com.tw)

# karma

[www.KarmaMedical.com](http://www.KarmaMedical.com)