


Ergo Aqua

Ergo Aqua is a two-in-one commode and shower chair designed to provide a premium lifestyle in everyday home care. It inherits the lightweight and ergonomic comfort of the Ergo Series.



The add-on EVA seat cushion ▲ is both comfortable and non-slip, ensuring a secure and enjoyable showering experience every day.



Premium EVA: low odor retention, durable, and easy to clean.



Flip-back armrests for easy access into the shower chair.



Low seat height to accommodate petite user. Adjustable seat height ranges from 400 to 475mm.

Fold and Stand ▶

With just a simple tap of your foot, Ergo Aqua effortlessly folds into an upright position for easy storage.



◀ Two-way Detachable Commode



Exclusive Splash-Proof Commode

The oval design provides a seamless fit under the seat cushion, eliminating gaps and effectively containing body waste.



Ergo Aqua Specifications

Overall L x W x H (mm)	580-605 x 520 x 800-875
Folded L x W x H (mm)	420 x 520 x 975
Seat Height (mm)	400-475
Max. User Weight (kg)	100
Total Mass (kg)	6.4

Remark : Karma reserves the right to modify information herein without further notice. " / " means optional ; " - " means adjustable ± 1 cm ± 0.5kg
Seat height and total mass specifications not included shower chair pad.
Rev.00 Aug. 2024



Scan to learn more about Ergo Aqua

KARMA MEDICAL PRODUCTS CO., LTD.

T +886-5-206 6688

F +886-5-206 7788

Email : globalsales@karma.com.tw

karma
www.KarmaMedical.com