

Ergo Lite

Reliable ultra-lightweight travel

The Ergo Lite is an ultra-lightweight wheelchair that weighs a mere 8.6 kg complete.

The compact design is ideal for transportation and storage. It is equipped with Karma's unique S-Ergo seating system. The Ergo Lite is an ideal choice for those who need a simple, lightweight wheelchair.



Ergo Lite



Ergo Lite 2



S-Ergo Seating Design (Shaped to fit you)

The revolutionary S-Ergo system is the first s-shaped ergonomic seating in the world, designed according to the body's curve, which is capable of simultaneously relieving pressure, stabilizing seating position, and reduces the tendency to slide forward.

Ultralight and Sturdy

We strive to make the caretaker's life easier. Weighing at a mere 8.6 kg, the ultralight weight lowers the risk of back injuries by reducing the physical effort required to lift it. Sturdy and safe, we believe it is the lightest crashed tested wheelchair available on the market.



Desk Armrests

The armrests are designed to allow the user to get closer to a desk-top.



Ergo Lite

KM-2501



Attendant-propelled wheel




Push & Pull Brake System

The dual function braking system allows carers to slow or stop the wheelchair.

Foldable Backrest and Frame

Allows for convenient storage and transportation. The foldable backrest and frame, makes storage and transport much easier.



Robust and Reliable

Ergo Lite is approved for transportation in accordance with ISO 7176-19 Crash Test. This means it is suitable to be used as a seat in a motor vehicle for transportation.

Karma Mobility Limited
S12668
20th Jan 2013





Ergo Lite 2

KM-2512

Improved functions

The Ergo Lite 2 has the basic functions of the Ergo Lite but is improved with 20" rear wheels and detachable footrests.



Self-propelled wheel



Detachable Footrests

Swing in and out detachable footrests add convenience when the user transfers in tight spaces.



Compact and Transportable

After easily removing the rear wheels and footrests, the transportable piece weighs only 9.8 kg. This makes storage and transport much easier.



Attendant-propelled wheel

Optional Accessories



Cup Holder



Travel bag



Pelvic Belt



Anti-tippers



Foldable Headrest



Net Bag



Backrest Bag



Quick Release Wheels
Only for Ergo Lite 2





Frame Color



Pearl Silver

Upholstery Color



Gray & Black

Specifications	Ergo Lite	Ergo Lite2
Seat Width* Seat Depth (inch)	16*17 / 18*17	16*17 / 18*17
Effective Seat Width (mm)	400 / 450	400 / 450
Effective Seat Depth (mm)	430	430
F/R wheels (inch)	F: 6 R: 14	F: 6 R: 14 / 20F / 20Q
Seat Surface Height at Front Edge (mm)	480	480
Backrest Height (mm)	420	420
Armrest Height (mm)	210	215
Calf Length (mm)	395-410	360-410
Max. User Weight (kg)	100	100
Total Mass (kg)	8.6 / 8.8	9.2-10.9
Overall L x W x H (mm)	930*560*905/ 930*610*905	min. 935*570*900 max. 930*645*900
Folded L x W x H (mm)	830*275*700/ 830*280*710	min. 635*280*700 max. 695*295*710

Remark : Karma reserves the right to modify information herein without further notice.
 " / " means optional ; " - " means adjustable ± 1 cm ± 0.5kg Rev.00 Jan.2021



KARMA MEDICAL PRODUCTS CO., LTD.

T +886-5-206 6688

F +886-5-206 7788

Email : globalsales@karma.com.tw

karma

www.KarmaMedical.com