

# Ergo Traveler

Ergo Traveler is an affordable rear-wheel drive power wheelchair featuring Karma's patented S-Ergo Seating System, which provides comfort and stability while promoting proper seating posture. The compact design makes Ergo Traveler ideal for confined spaces and short-distance travel.



## S-Ergo Seating System

The unique seating system can release pressure, stabilize posture, and reduce sliding.



## Indoor Maneuverability

Smooth and easy to navigate, this chair moves effortlessly through doorways, onto lifts, or around confined spaces.



## Foldable and Travel Friendly

Disassemble and reassemble with no tools required, making travel a breeze!

## Adjustable and Transfer Friendly

- Height adjustable and flip-up armrests: extra comfort and effortless transfer in and out of the chair.
- Detachable swing-away footrests: can turn outward or detach for effortless transfer in and out of the chair.

## OPTIONS

- Swing-away controller bracket
- Quick release attendant control
- Attendant brake

Nominal seat width (mm)	Nominal seat depth (mm)	Effective seat width (mm)	Effective seat depth (mm)	Backrest height w/Cushion (mm)	Seat height w/o Cushion (mm)	Seat height w/Cushion (mm)	Power seat lift (mm)	Power Tilt, Seat Plane (X°)	Backrest recline, seat to back angle (X°)	Legrest angle, leg to seat (X°)	F/M/R Wheel (in)	Overall length with legrest (mm)	Overall width (mm)
460 (18)	460	450	475	500	-	540	-	-	88-102	90	8/-/12	1030-1075	620

Max. User weight (kg)	Drive range (km)	Max. speed (km)	Obstacle climbing (mm)	Min. turning radius w/o user (mm)	Ability to climb rated slope (X°)	Battery (Ah)	Electrical system
120	14.9(22 Ah) 18.7(36 Ah)	6	50	755	6	22*2 36*2	VR2

- The actual product specifications may differ from the data in the table by ±10 mm ±0.5 kg.  
 - The actual product specifications may vary according to different configurations.  
 - Karma reserves the right to modify information herein without further notice.

