

Falcon

Falcon features a foldable frame and a removable battery pack, giving users full-size performance in a lightweight and compact package.



Compact and Maneuverable

A compact medium-size wheelchair, with an overall width of 590 mm, makes moving around tight corners and through narrow doorways a breeze.



Height Adjustable Armrests

The armrest height can be adjusted to four different levels according to the user's body requirements.



Easy Transfer

Flip-up armrests and detachable swing-away footrests allow for both lateral and forward transfers.



Tool-free Quick Folding

The lightweight, foldable frame and the quick-release battery design make this wheelchair easy to transport and store.



OPTIONS

- Swing-away controller bracket
- Quick release attendant control
- Foldable headrest
- Elevating legrest
- Handle bar

Nominal seat width (mm)	Nominal seat depth (mm)	Effective seat width (mm)	Effective seat depth (mm)	Backrest height w/Cushion (mm)	Seat height w/o Cushion (mm)	Seat height w/Cushion (mm)	Power seat lift (mm)	Power Tilt, Seat Plane (X°)	Backrest recline, seat to back angle (X°)	Legrest angle, leg to seat (X°)	F/M/R Wheel (in)	Overall length with legrest (mm)	Overall width (mm)
410/460/515	430/480/480	430/480/530	470/515/515	450-525	-	525	-	-	95	110	8/-/12	985-1115	590/640/690

Max. User weight (kg)	Drive range (km)	Max. speed (km)	Obstacle climbing (mm)	Min. turning radius w/o user (mm)	Ability to climb rated slope (X°)	Battery (Ah)	Electrical system
115/135 (pneumatic tires)	32	6	50	650/825/1000	12	50*2	VR2

- The actual product specifications may differ from the data in the table by ±10 mm ±0.5 kg.
 - The actual product specifications may vary according to different configurations.
 - Karma reserves the right to modify information herein without further notice.

