

Saber

Karma's Saber series has provided users with unparalleled outdoor performance for over 20 years. World-renowned for its speed, handling, and reliability, Saber is the ideal power wheelchair for users with good upper body strength and a taste for adventure.



Easy Transfer

Locking flip-up armrests offer support and security while driving and convenience while transferring.



Comfort and Support

Optional power tilt and recline functions allow users to redistribute and relieve pressure, while the standard handlebar provides extra stability.

The Legend of Mr. Choi

From 2006-2007, South Korean Chang-Hyun Choi, who has cerebral palsy, completed a 28,000 km journey across 35 countries with his mouth-controlled Saber.

Source: Guinness World Records



Born for the Outdoors



Illuminate your path with Saber's LED lighting system.



Reflective stripes add safety with style, while the push bar allows attendants to comfortably maneuver Saber in freewheel mode.



Independent Suspension Response (ISR) acts to absorb bumps, maximizing comfort and stability.

Water-resistant accessories and all-weather tires prepare Saber for any journey!



Saber

FRAME COLORS

-  Ruby Red
-  Navy Blue
-  Dark Green

Model	KP-45.5	KP-45.5T	KP-45.5TR
Driving Range (km)	41 (50 Ah) / 64 (80 Ah)	41 (50 Ah)	41 (50 Ah)
Manual Recline	0-151	88-110	
Power Tilt-in-Space		6-45	7-45
Power Recline			90-119

Product Name	Nominal seat width (mm)	Nominal seat depth (mm)	Effective seat width (mm)	Effective seat depth (mm)	Backrest height w/Cushion (mm)	Seat height w/o Cushion (mm)	Seat height w/Cushion (mm)	Power seat lift (mm)	Power Tilt, Seat Plane (X°)	Backrest recline, seat to back angle (X°)	Legrest angle, leg to seat (X°)	F/M/R Wheel (in)	Overall length with legrest (mm)
Saber	460/480	425/485	470-500/490-570	510/505	535/560	-	575/550	-	-	0-147/0-151	107-179	9/-/16	1100-1670/1085-1540
Saber T	480	480	480-500	510	550	-	580	-	6-45	88-110	102-179	9/-/16	1120-1430
Saber TR	480	480	480-500	510	550	-	580	-	7-45	90-119	102-179	9/-/16	1120-1600

Overall width (mm)	Max. User weight (kg)	Drive range (km)	Max. speed (km)	Obstacle climbing (mm)	Min. turning radius w/o user (mm)	Ability to climb rated slope (X°)	Battery (Ah)	Electrical system
670-690/705-730	135	40 (50Ah)/64 (80Ah)	13	50	945	12	50*2/80*2	VR2/Linx
670-690	135	40	13	50	945	12	50*2/80*2	VR2/Linx
670-690	135	40	13	50	945	12	50*2/80*2	VR2/Linx

- The actual product specifications may differ from the data in the table by ±10 mm ±0.5 kg.
 - The actual product specifications may vary according to different configurations.
 - Karma reserves the right to modify information herein without further notice.

